



## 5 reasons to choose Warwick SASCo


Recognised as the market leader in the design and  
manufacture of medical polyware

Wide range of high quality re-usable and single-use polypropylene  
products available

Unrivalled reputation for providing quality products, delivery and  
service to healthcare professionals worldwide

Re-usable products are guaranteed for a minimum of 1000 wash  
cycles – delivering excellent value for money

Certification since 1997 demonstrates our continued commitment  
to maintaining industry standards. Our quality management system  
complies to ISO13485 and ISO9001



 **MADE IN BRITAIN**



## Warwick SASCo

For over 30 years Warwick SASCo has been part of evolving supply of decontamination and nursing products worldwide to hospitals, healthcare and patient nursing.

Recognised as a market leader, we have a reputation for high quality products designed for hospital decontamination and sterile services, operating and surgical applications, patient care and home nursing - where infection management and good user application is paramount.

Warwick SASCo's low cost reusable and single use product ranges are chosen by healthcare professionals and end users worldwide.



Our growing family owned business is proud to be represented by leading distributor partners in over 60 countries.

All our products are manufactured within the UK with Quality Assured Suppliers. Our customer support and distribution centre is in Warwick, England.

The team at Warwick SASCo is immensely proud of its reputation for delivering high quality products and a professional friendly service to you - our customer. You will find our experienced team are knowledgeable, caring and quick to respond to your requirements.

A handwritten signature in blue ink that reads "Darby Booth". The signature is written in a cursive style and is positioned above the printed name.

Darby Booth,  
Managing Director – Warwick SASCo



# Products Overview

Our wide range of reusable and disposable medical polyware products are chosen by healthcare professionals worldwide for their high quality and low cost.

 <h2>Reusable Medical Polyware</h2> <p>WARWICK RANGE</p>	<p>The Warwick Range for use in Decontamination Services, Operating Rooms, Surgeries, Clinics and healthcare services. Products are designed for repeated washing, disinfection and sterilization.</p> <p>The Warwick Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).</p> <p>Products comply with BS 5452 for Hospital Hollowware where applicable. Material used is polypropylene unless otherwise stated.</p>
 <h2>Single Use Medical Polyware</h2> <p>STERILITE RANGE</p>	<p>Sterilizable lightweight polyware for single-use and disposable sets. SteriLite Range is designed for Custom Pack Manufacturers and Decontamination Services.</p> <p>The Sterilite Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).</p> <p>Material used is polypropylene unless otherwise stated.</p>
 <h2>Nursing &amp; Care Daily Living Aids &amp; Accessories</h2> <p>SPA RANGE</p>	<p>Spa Range of personal living aids, patient care, toileting and bedside products.</p> <p>Products comply with BS 2588 for Portable Sanitary Products and BS 3215 for Reusable Urinals where applicable.</p> <p>Material used is polypropylene unless otherwise stated.</p>



## How to Order

To place an order, request a quotation or for general enquiries:

CALL us on +44 (0)1926 422 427

FAX +44 (0)1926 315 405


EMAIL [enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

We work with a wide range of industry sectors including:



**STERILE SERVICES & DECONTAMINATION**

Operating Room,  
Instrument Packing




**PERSONAL USE & HOME CARE**

Daily Living Aids,  
Community Nurse,  
Incontinent Care



**HOSPITALS**

Surgery, Clinics, Wards,  
Veterinary, Dentistry



**NURSING HOMES**

Patient Care,  
Rehabilitation



**CUSTOM PACK MANUFACTURERS**

Single Use Packs,  
Contract Instrument  
Processing



To find out more about our service, product ranges, product decontamination cycle and industry news:

visit [www.sasco.co.uk](http://www.sasco.co.uk)

**WARWICK - reusable range**

Compartment Trays	page	6
DIN Style Baskets	page	7-8
Disinfectant & Soaking Containers	page	6
Gallipots	page	9-10
Instrument Trays	page	7-8
Instrument Tray Lids	page	9
Jugs	page	10
Kidney Dishes	page	11
Lotion & Sponge Bowls	page	9
Medical Boxes, Instrument Boxes & Storage Containers	page	12
Scope Containers	page	13
Silicone Protection	page	13
Theatre Bowls	page	11
Tray Tags	page	14
Wash Bowls	page	11

**STERILITE - single use range**

Compartment Trays	page	16
Gallipots	page	17
Jugs	page	18
Kidney Dishes	page	18
Packing Trays	page	16
Quivers	page	18
Sponge Bowls	page	17
Theatre Bowls	page	18

**SPA - care and living aid range**

Bedpans	page	20
Bedpan Holders	page	21
Bidets	page	20-21
Commode Pans	page	20-21
Denture Cups	page	21
Dispensing and Dressing Trays	page	23
Drinking Beakers and Cups	page	22
Drinking Cup Lids	page	22
Jugs	page	24
Jug Sets	page	24
Kidney Dishes	page	23
Medicine Measures	page	24
Urinal Bottles	page	25
Urinal Pans	page	25
Vomit Bowls	page	23

# WARWICK RANGE

## Reusable Medical Polyware

The Warwick Range is a wide range of sterilizable long-life plastic polyware products.

Designed as an alternative to stainless steel products, they are ideal for use in Decontamination Services, Operating Rooms, Clinics, Surgeries and other surgical or healthcare services.

Products are designed for repeated washing, disinfection and sterilization. The Warwick Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).

Products comply with BS 5452 for Hospital Hollowware where applicable.

Material used is polypropylene unless otherwise stated.

## Warwick Benefits

- Low cost.
- Robust and hard wearing - guaranteed for 1000 cycles.
- Easy to clean - blood protein is easily removed from the surface, leaving no residue.
- 50% lighter than stainless steel, assisting users with storage, lifting and handling.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).
- Can be used as standalone items or as a support component in surgical sets.
- Widely used in Sterile Services for packing instrument sets and procedure kits for the autoclave.
- Disinfection sets can be used for manual sterilization in Sterile Services, Clinics, etc.
- Used non-sterile in hospital wards and clinics.
- Smooth corners and surfaces, reducing risk of paper rips and the need for re-processing.
- Autoclave sterilization temperature is 134.6°C + 2°C (generally referred to as 135°C).



# Compartment Trays, Disinfectant and Soaking Containers



Compartment Tray -  
3 compartments

**Code: CIT1814**

Dimensions: 181x134x54mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 38

R W S E HP LF



Compartment Tray -  
3 compartments

**Code: CIT2013**

Dimensions: 205x130x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 120

R W S E HP LF



Compartment Tray -  
4 compartments

**CIT2718**

Dimensions: 270x180x41mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF



Compartment Tray -  
4 compartments

**Code: LCT2215**

Dimensions: 220x150x24mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 160

R W S E HP LF



Disinfection Set -  
Soaking container, strainer  
and hinged lid

**Code: DSSW1000**

Dimensions: 260x130x65mm  
Volume: 1litre  
Material: Polypropylene  
Colour: White  
Pack Size: 10

R W S E HP LF



Disinfection Set -  
Soaking container, strainer  
and lid

**Code: DSSW2711**

Dimensions: 257x180x128mm  
Volume: 3litre  
Material: Polypropylene  
Colour: White  
Pack Size: 5

R W S E HP LF



Disinfection Set -  
Soaking container, strainer  
and lid

**Code: DSS2711**

Dimensions: 257x180x128mm  
Volume: 3litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 5

R W S E HP LF



Disinfection Set -  
Soaking container, strainer  
and vented lid

**Code: DSSV2711**

Dimensions: 257x180x128mm  
Volume: 3litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 5

R W S E HP LF



Disinfection Set -  
Soaking container, strainer,  
drip tray and lid

**Code: DSS4215**

Dimensions: 426x156x138mm  
Volume: 4litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 1

R W S E HP LF



Disinfection Set -  
Soaking container, strainer,  
drip tray and lid

**Code: DSS8515**

Dimensions: 890x200x145mm  
Volume: 12litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 1

R W S E HP LF



Disinfection Set - 2 x soaking  
containers, 1 x strainer, 2 x lids  
and 1 x drip tray

**Code: DSS8515TP**

Dimensions: 890x200x145mm  
Volume: 12litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 1

R W S E HP LF

## How to order

To place an order,  
request a quotation or for  
general enquiries:

CALL us on  
**+44 (0)1926 422 427**

FAX

**+44 (0)1926 315 405**

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

# DIN Style Baskets and Instrument Trays



Instrument Tray -  
Solid base

**Code: IT1813**

Dimensions: 181x134x54mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 40



Instrument Tray -  
Solid base

**Code: IT1913**

Dimensions: 195x135x23mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 115



Instrument Tray -  
Solid base

**Code: T2015**

Dimensions: 200x150x51mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60



Instrument Tray -  
Solid ribbed base

**Code: IT2015**

Dimensions: 200x150x51mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60



Instrument Tray -  
Solid base

**Code: IT2715**

Dimensions: 270x150x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60



Instrument Tray -  
Solid base

**Code: T3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25



## How to order

To place an order,  
request a quotation or for  
general enquiries:

CALL us on  
**+44 (0)1926 422 427**

FAX  
**+44 (0)1926 315 405**

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)



Instrument Tray -  
Solid ribbed base

**Code: IT3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25



Instrument Tray -  
Solid ribbed base

**Code: IT4030**

Dimensions: 424x305x75mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 12



**R** Reusable

**W** Washer Decontamination

**S** Suitable for Steam Sterilization

**E** Suitable for Ethylene Oxide Sterilization

**HP** Suitable for Hydrogen Peroxide Sterilization

**LF** Latex Free



# DIN Style Baskets and Instrument Trays



Instrument Tray -  
Mesh base and sides

**Code: MIT1515**

Dimensions: 200x180x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 50

R W S E HP LF



Instrument Tray -  
Mesh base

**Code: MIT2015**

Dimensions: 200x150x51mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF



Instrument Tray -  
Mesh base and sides

**Code: MIT2709**

Dimensions: 275x110x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF



Instrument Tray -  
Mesh base

**Code: MIT2715**

Dimensions: 270x150x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF



Instrument Tray -  
Mesh base and sides

**Code: MIT2817**

Dimensions: 280x170x55mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF



Instrument Tray -  
Mesh base and sides

**Code: MIT3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 24

R W S E HP LF



Instrument Tray -  
Mesh base

**Code: MIT4030**

Dimensions: 424x305x75mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 12

R W S E HP LF



Instrument Tray -  
Perforated base

**Code: PIT3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF



Instrument Tray -  
Perforated base

**Code: PIT4030**

Dimensions: 424x305x75mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 12

R W S E HP LF



DIN Style Tray -  
Mesh base and sides

**Code: DINT2525**

Dimensions: 255x255x50mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 20

R W S E HP LF



DIN Style Tray -  
Mesh base and sides

**Code: DINT4825**

Dimensions: 480x255x50mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

# Instrument Tray Lids



Instrument Tray Lid -  
For trays IT1813/CIT1814

**Code: ITL181**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10



Instrument Tray Lid -  
For trays IT2015/MIT2015/  
T2015

**Code: ITL201**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10



Instrument Tray Lid -  
For trays IT2715/MIT2715

**Code: ITL271**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10



Instrument Tray Lid -  
For tray CIT2718

**Code: ITL281**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10



Instrument Tray Lid -  
For tray MIT2817

**Code: ITL287**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10



Instrument Tray Lid -  
For trays IT3025/MIT3025/  
PIT3025/T3025

**Code: ITL302**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10



## How to order

To place an order, request a quotation or for  
general enquiries:

CALL us on

+44 (0)1926 422 427

FAX

+44 (0)1926 315 405

EMAIL

[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

# Gallipots, Lotions and Sponge Bowls



Gallipot

**Code: GP40**

Dimensions: 40x30mm  
Volume: 35ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 100



Gallipot

**Code: GP60**

Dimensions: 60x40mm  
Volume: 95ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 300



Gallipot

**Code: GP80**

Dimensions: 80x60mm  
Volume: 280ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 80



Gallipot -  
Flat base

**Code: GPLB60**

Dimensions: 60x40mm  
Volume: 50ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 400



# Gallipots, Lotion and Sponge Bowls



Gallipot - Flat base

Code: GPLB80

Dimensions: 80x60mm  
Volume: 200ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 200

R W S E HP LF



Lotion Bowl - Flat base

Code: GPLB100

Dimensions: 100x45mm  
Volume: 300ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 135

R W S E HP LF



Lotion Bowl - Flat base

Code: GPLB150

Dimensions: 150x70mm  
Volume: 500ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF



Sponge Bowl - Flat base

Code: SB140

Dimensions: 140x65mm  
Volume: 300ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 160

R W S E HP LF



Lotion Bowl - Graduated

Code: LB100

Dimensions: 100x45mm  
Volume: 280ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 200

R W S E HP LF



Lotion Bowl

Code: LB150

Dimensions: 150x70mm  
Volume: 900ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF



Lotion Bowl

Code: LB200

Dimensions: 200x90mm  
Volume: 2000ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 30

R W S E HP LF



Lotion Bowl

Code: LB250

Dimensions: 250x105mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 20

R W S E HP LF

# Jugs



Jug

Code: JG500

Dimensions: N/A  
Volume: 500ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 50

R W S E HP LF



Jug

Code: JG1000

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 45

R W S E HP LF



Jug

Code: JG2000

Dimensions: N/A  
Volume: 2000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 20

R W S E HP LF



# Kidney Dishes, Theatre Bowls and Wash Bowls



Kidney Dish

**Code: KD150**

Dimensions: 150x35mm  
Volume: 300ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 210



Kidney Dish

**Code: KD200**

Dimensions: 200x45mm  
Volume: 500ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 150



Kidney Dish

**Code: KD250**

Dimensions: 250x55mm  
Volume: 750ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60



Kidney Dish

**Code: KDY250**

Dimensions: 250x55mm  
Volume: 750ml  
Material: Polypropylene  
Colour: Yellow  
Pack Size: 60



Kidney Dish

**Code: KD300**

Dimensions: 300x65mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25



Kidney Dish Tray -  
Shallow kidney shaped tray

**Code: KDT185**

Dimensions: 185x92mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 70



Washbowl -  
Graduated

**Code: WB3112**

Dimensions: 310x105mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25



Wash Bowl

**Code: WBN3112**

Dimensions: 310x105mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 25



Theatre Bowl

**Code: WB3513**

Dimensions: 350x 20mm  
Volume: 5000ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25



Wash Bowl

**Code: WBN3513**

Dimensions: 350x 20mm  
Volume: 5000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 25



**R** Reusable

**W** Washer Decontamination

**S** Suitable for Steam Sterilization

**E** Suitable for Ethylene Oxide Sterilization

**HP** Suitable for Hydrogen Peroxide Sterilization

**LF** Latex Free

# Medical Boxes, Instrument Boxes & Storage Containers



Instrument Box -  
Storage box with lid

Code: IBX2718

Dimensions: 257x180x128mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 20

R W S E HP LF



Instrument box -  
Storage box with vented lid

Code: IBXV2718

Dimensions: 257x180x128mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 20

R W S E HP LF



Instrument Box -  
Storage box with lid

Code: IBXN290

Dimensions: 300x224x140mm  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 12

R W S E HP LF

## How to order

To place an order,  
request a quotation or for  
general enquiries:

CALL us on  
+44 (0)1926 422 427

FAX  
+44 (0)1926 315 405

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)



Transportation Box -  
With hinged lid

Code: MB2318G

Dimensions: 230x170x200mm  
Volume: N/A  
Material: ABS  
Colour: Green  
Pack Size: 10

R W S E HP LF



Transportation Box -  
With hinged lid

Code: MB2318

Dimensions: 230x170x200mm  
Volume: N/A  
Material: ABS  
Colour: Orange  
Pack Size: 10

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

# Scope Containers



Scope Container -  
With silicone dividers and lid

**Code: SMT5510**

Dimensions: 660x110x66mm  
Volume: N/A  
Material: Polypropylene  
Colour: Lilac  
Pack Size: 1

R W S E HP LF



Silicone Fixer -  
For SMT5510

**Code: SF9225**

Dimensions: N/A  
Volume: N/A  
Material: Silicone  
Colour: Blue  
Pack Size: 10

R W S E HP LF



Silicone Plug -  
For SMT5510

**Code: SG5510**

Dimensions: N/A  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF

# Silicone Protection



Silicone Pad -  
For trays IT/MIT2715 and  
small metal DIN tray  
**Code: SP2313**

Dimensions: 230x132mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF



Silicone Pad -  
For trays IT/MIT/PIT/T3025

**Code: SP2823**

Dimensions: 280x230mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF




Silicone Pad -  
Fits 1/2 tray IT/MIT/PIT4030

**Code: SP2719**

Dimensions: 270x190mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF



Silicone Pad -  
Fits 1/2 tray DINT4825

**Code: SP4622**

Dimensions: 230x220mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF



Silicone Pad -  
For medium metal DIN tray

**Code: SP2324**

Dimensions: 235x245mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF




Silicone Pad -  
For large metal DIN tray

**Code: SP4724**

Dimensions: 470x245mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF



Silicone Roll -  
1 metre length with 10  
divisible sections

**Code: SR10**

Dimensions: 100x245mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 1

R W S E HP LF

# Tray Tags



Tray Tag

**Code: TAG-010**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: White  
Pack Size: 100



Tray Tag

**Code: TAG-015**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Yellow  
Pack Size: 100



Tray Tag

**Code: TAG-020**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Blue  
Pack Size: 100



Tray Tag

**Code: TAG-022**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Light Green  
Pack Size: 100



Tray Tag

**Code: TAG-025**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Orange  
Pack Size: 100



Tray Tag

**Code: TAG-030**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Red  
Pack Size: 100



Tray Tag

**Code: TAG-050**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Green  
Pack Size: 100



Tray Tag

**Code: TAG-055**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Grey  
Pack Size: 100



Tray Tag

**Code: TAG-060**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Lavender  
Pack Size: 100



Tray Tag

**Code: TAG-065**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Tan  
Pack Size: 100



Tray Tag

**Code: TAG-070**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Black  
Pack Size: 100



Stainless Steel Attachment Ring

**Code: TR30**

Dimensions: 30mm  
Volume: N/A  
Material: Steel  
Colour: Silver  
Pack Size: 100



**R** Reusable

**W** Washer Decontamination

**S** Suitable for Steam Sterilization

**E** Suitable for Ethylene Oxide Sterilization

**HP** Suitable for Hydrogen Peroxide Sterilization

**LF** Latex Free

# STERILITE RANGE

## Single Use Medical Polyware

The SteriLite Range is a wide range of single use plastic polyware products.

Designed to speed up packing and sterilizing procedures. SteriLite is ideal for use in Sterile Pack Manufacturing, Decontamination Services, Operating Rooms and Clinics due to its low cost, product design and performance in hospitals.

The Sterilite Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).

Material used is polypropylene unless otherwise stated.

## SteriLite Benefits

- Low cost and disposable.
- No need to invest time and money into retrieval and decontamination after use.
- Clean and ready for sterilization.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).
- Can be used as stand-alone items or as a support (supplementary) component in operating rooms and surgical procedures.
- Widely used in Sterile Services for packing instrument sets and procedure kits.
- Used non-sterile in hospital wards and clinics.





# Compartment Trays and Packing Trays



Packing Tray -  
4 compartments

**Code: LCT2215S**

Dimensions: 220x150x24mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 160




Packing Tray

**Code: ST4330**

Dimensions: 430x330x64mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 63



 Single Use

 Washer Decontamination

 Suitable for Steam Sterilization

 Suitable for Ethylene Oxide Sterilization

 Suitable for Hydrogen Peroxide Sterilization

 Latex Free

# Gallipots and Sponge Bowls



Gallipot

**Code: GPS60**

Dimensions: 62x30mm  
 Volume: 60ml  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 2000

- Single Use
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free



Gallipot

**Code: GPS120**

Dimensions: 85x30mm  
 Volume: 120ml  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 1000

- Single Use
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free



Gallipot

**Code: GP100**

Dimensions: 100x53mm  
 Volume: 250ml  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 750

- Single Use
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free



Gallipot

**Code: GP130**

Dimensions: 130x66mm  
 Volume: 500ml  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 735

- Single Use
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free



Sponge Bowl

**Code: SB170**

Dimensions: 170x84mm  
 Volume: 1000ml  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 340

- Single Use
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free



Sponge Bowl

**Code: SB200**

Dimensions: 210x90mm  
 Volume: 2000ml  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 220

- Single Use
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free

- Single Use
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free

## Jugs and Quivers



Jug

Code: J1000

Dimensions: N/A  
 Volume: 1000ml  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 100



Quiver

Code: QV30

Dimensions: 190x65mm  
 Volume: N/A  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 75



## Kidney Dishes and Theatre Bowls



Kidney Dish

Code: KB250

Dimensions: 252x100x48mm  
 Volume: 750ml  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 450



Kidney Dish

Code: KB300

Dimensions: 300x140x57mm  
 Volume: 1500ml  
 Material: Polypropylene  
 Colour: Natural  
 Pack Size: 135



Theatre Bowl

Code: TB350

Dimensions: 352x116mm  
 Volume: 5000ml  
 Material: Polypropylene  
 Colour: Blue  
 Pack Size: 50



### How to order

To place an order,  
 request a quotation or for  
 general enquiries:

CALL us on  
 +44 (0)1926 422 427  
 FAX  
 +44 (0)1926 315 405

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

- Single Use
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free

# SPA RANGE

## Nursing & Care, Daily Living Aids and Accessories

Spa Range is a wide range of re-usable personal living aids, patient care, toileting and bedside products.

The Spa Range is designed for use in Hospital Wards, Rehabilitation, Nursing Homes, Community Care and for Personal Use in the Home.

Products comply with BS 2588 for Portable Sanitary Products and BS 3215 for Reusable Urinals where applicable.

Material used is polypropylene unless otherwise stated.

### Spa Benefits

- Toileting products that preserve patient dignity and independence by empowering patients to toilet themselves with little or no assistance.
- Urinal Bottles specifically designed to the male and female anatomy for a comfortable and safe fit during use.
- Product design reduces risk of cross infection, unpleasant odours and sights.
- Easy to use, clean and disinfect, meaning both patient and carer can be confident that products are safe to re-use.
- Beakers and Drinking Cups are wide, easy to grip, non-scratch, drop-proof and have a nondrip rim.
- Handled beakers are ideal for patients with unsteady or infirm hands.
- Medicine Measures are sturdy, ideal for long-term use and regularly checked by Trading Standards for measurement accuracy and compliance.
- All Spa products can be cleaned in automatic washer disinfectors. Also suitable for hand washing using detergents and disinfectants.
- Urinal Pan SPN and Bedpans FPN, VPN, DPNF, BPN and CMP can be autoclaved.



# Bedpans, Commode Pans and Bidets



Bedpan -  
Adult

**Code: BPN**

Dimensions: 370x310x90mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 30



Slipper Pan -  
Bedpan with handle and lid

**Code: FPN**

Dimensions: 449x309x99mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 12



Slipper Pan -  
Bedpan with handle and lid

**Code: FPNW**

Dimensions: 449x309x99mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: White  
Pack Size: 12



Bedpan -  
Hospital style bedpan

**Code: VP**

Dimensions: 355x308x108mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 18



Bedpan -  
Hospital bedpan with lid

**Code: VPN**

Dimensions: 355x308x108mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 15



Bedpan -  
Hospital bedpan with lid

**Code: VPNW**

Dimensions: 355x308x108mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: White  
Pack Size: 15



Bedpan -  
European bedpan with handle  
and lid

**Code: DPNF**

Dimensions: 420x315x100mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 15



## How to order

To place an order,  
request a quotation or for  
general enquiries:

CALL us on  
**+44 (0)1926 422 427**  
FAX  
**+44 (0)1926 315 405**

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)



Slipper Pan -  
Female urinal

**Code: SPN**

Dimensions: 335x245mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 25



Commode Pan / Bidet

**Code: CMP**

Dimensions: 359x338x130mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 20



Commode Pan / Bidet

**Code: CMPW**

Dimensions: 359x338x130mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: White  
Pack Size: 20



Commode Pan -  
Commode pan with lid

**Code: CBL**

Dimensions: 259x259x144mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: White  
Pack Size: 25



**R** Reusable

**W** Washer Decontamination

**S** Suitable for Steam Sterilization

**E** Suitable for Ethylene Oxide Sterilization

**HP** Suitable for Hydrogen Peroxide Sterilization

**LF** Latex Free

## Bedpan Holders, Commode Pans and Bidets



Bedpan Support -  
For pulp disposable bedpan

**Code: PBPS**

Dimensions: 370x310x90mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Cream  
Pack Size: 30



Bedpan Support -  
For pulp disposable midi  
slipper pan

**Code: PSPS**

Dimensions: 470x290x75mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Cream  
Pack Size: 30



## Denture Cups

Denture Cup -  
With strainer and lid

**Code: DNT**

Dimensions: 90x80x60mm  
Volume: N/A  
Material: Polypropylene  
Colour: Green  
Pack Size: 75



### How to order

To place an order,  
request a quotation or for  
general enquiries:

CALL us on  
+44 (0)1926 422 427

FAX  
+44 (0)1926 315 405

EMAIL  
enquiries@sasco.co.uk

# Drinking Beakers, Cups and Lids



Beaker -  
Feeder/drinking cup

Code: BFC250

Dimensions: N/A  
Volume: 250ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 300

R W S E HP LF



Beaker -  
Feeder/drinking cup with  
handles

Code: BFH250

Dimensions: N/A  
Volume: 250ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 140

R W S E HP LF



Lid -  
Narrow Spout for BFC250/  
BFH250

Code: BNS

Dimensions: 4x4mm spout  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 300

R W S E HP LF



Lid -  
Medium Spout for BFC250/  
BFH250

Code: BMS

Dimensions: 7x4mm spout  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 300

R W S E HP LF



Lid -  
Wide spout for BFC250  
& BFH250

Code: BWS

Dimensions: 13x8mm spout  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 300

R W S E HP LF



Lid -  
Rim feeder for beakers BFC250  
& BFH250

Code: BRF

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 180

R W S E HP LF



Drinking Cup -  
With handles, spout lid and  
feeder lid

Code: GPMC6250

Dimensions: N/A  
Volume: 250ml  
Material: Polypropylene  
Colour: Cream  
Pack Size: 90

R W S E HP LF

## How to order

To place an order,  
request a quotation or for  
general enquiries:

CALL us on

+44 (0)1926 422 427

FAX

+44 (0)1926 315 405

EMAIL

enquiries@sasco.co.uk



Beaker Set -  
Feeder/drinking cup with han-  
dle, spout lid and feeder lid

Code: HDMS

Dimensions: N/A  
Volume: 300ml  
Material: SAN  
Colour: Clear  
Pack Size: 40

R W S E HP LF



Teapot Style Beaker

Code: TFB

Dimensions: N/A  
Volume: 250ml  
Material: Polypropylene  
Colour: White  
Pack Size: 45

R W S E HP LF



Clear Glass Style Beaker

Code: TBC

Dimensions: N/A  
Volume: 280ml  
Material: SAN  
Colour: Clear  
Pack Size: 120

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

# Dispensing and Dressing Trays, Kidney Dishes and Vomit Bowls



Dispensing/Dressing Tray

**Code: ITC2715**

Dimensions: 270x150x30mm  
 Volume: N/A  
 Material: Polypropylene  
 Colour: Green  
 Pack Size: 60

**R** **W** **S** **E** **HP** **LF**



Kidney Dish

**Code: KDC250**

Dimensions: 250x130mm  
 Volume: 750ml  
 Material: Polypropylene  
 Colour: Green  
 Pack Size: 60

**R** **W** **S** **E** **HP** **LF**



Vomit Bowl -  
With handle

**Code: VB300**

Dimensions: 300x144mm  
 Volume: 1500ml  
 Material: Polypropylene  
 Colour: Green  
 Pack Size: 25

**R** **W** **S** **E** **HP** **LF**

**R** Reusable

**W** Washer Decontamination

**S** Suitable for Steam Sterilization

**E** Suitable for Ethylene Oxide Sterilization

**HP** Suitable for Hydrogen Peroxide Sterilization

**LF** Latex Free



## Jugs and Jug Sets



Jug with lid

**Code: JGL1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 40



Jug

**Code: JG500**

Dimensions: N/A  
Volume: 500ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 50



Jug

**Code: JG1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 45



Jug Set -  
Includes tray, carafe with lid  
and beaker  
**Code: LST**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: N/A  
Pack Size: 20



### How to order

To place an order,  
request a quotation or for  
general enquiries:

CALL us on  
**+44 (0)1926 422 427**  
FAX  
**+44 (0)1926 315 405**

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

## Medicine Measures



Medicine Measure -  
Graduated

**Code: MM50**

Dimensions: N/A  
Volume: 50ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 1000



Medicine Measure -  
Graduated

**Code: MM60**

Dimensions: N/A  
Volume: 60ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 1000



Medicine Measure -  
Graduated black lines

**Code: MM60BL**

Dimensions: N/A  
Volume: 60ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 1000



# Urinal Bottles and Urinal Pans



Urinal Bottle -  
Female urinal with handle

**Code: UF1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polyethylene  
Colour: Natural  
Pack Size: 30



Urinal Pan -  
Female Urinal

**Code: PUF1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polyethylene  
Colour: White  
Pack Size: 20



Urinal Bottle -  
Male urinal with handle

**Code: UMA1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polyethylene  
Colour: Natural  
Pack Size: 40



Urinal Bottle -  
Male urinal with handle,  
hospital style

**Code: UBS1000**

Dimensions: N/A  
Volume: 1200ml  
Material: Polyethylene  
Colour: Natural  
Pack Size: 45



Urinal Bottle -  
Male urinal with handle

**Code: UML1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polyethylene  
Colour: Natural  
Pack Size: 50



Urinal Bottle -  
Male urinal, hospital style

**Code: UM1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 25



Urinal Bottle -  
Male urinal, hospital style

**Code: UM650**

Dimensions: N/A  
Volume: 650ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 30



Non Leak Valve -  
For use with UMA1000 &  
UML1000

**Code: UMT**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: White  
Pack Size: 25



Urinal Bottle with UMT -  
Male urinal bottle with non  
leak valve

**Code: USP1200**

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene and  
Polyethylene  
Colour: N/A  
Pack Size: 10



Urinal Bottle Holder -  
For UBS1000, UML1000 and  
UM1000

**Code: UBH2002**

Dimensions: N/A  
Volume: 1200ml  
Material: Plastic Coated Wire  
Colour: White  
Pack Size: 10



## How to order

To place an order, request a quotation or for  
general enquiries:

CALL us on

+44 (0)1926 422 427

FAX

+44 (0)1926 315 405

EMAIL

[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

**R** Reusable





**W** Washer Decontamination

**S** Suitable for Steam Sterilization

**E** Suitable for Ethylene Oxide Sterilization

**HP** Suitable for Hydrogen Peroxide Sterilization

**LF** Latex Free



Vertrieb in der Schweiz:  
Distributeur pour la Suisse:



Novimed AG  
Heimstrasse 46 | 8953 Dietikon  
Telefon +41 44 743 40 60  
Telefax +41 44 743 40 65  
info@novimed.ch | novimed.ch

Warwick House, Heathcote Way, Warwick, CV34 6TE  
T: +00 44 (0) 1926 422 427  
F: +00 44 (0) 1926 315 405  
E: enquiries@sasco.co.uk

[www.sasco.co.uk](http://www.sasco.co.uk)